

27641 QA

8 Ch Foot Ankle Coil Check 2020-04-29 15:28:10

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_8Ch_FootAnkle_Sag_sn.edb.

Measured Series Loid #1:

Local 2.0.145268191

Measured Series Loid #2:

Local 2.0.145262920

Ordered Series Loid #1:

Local 2.0.145268191

Series Loid for Evaluation 1:

Local 2.0.145268191

Number of images in series: 1

Evaluate image #1

Protocol 005_8Ch_FootAnkle_Sag_sn

No Evaluation for FA1_FA2_FA3_FA4_FA5_FA6_FA7_FA8 specified in ServiceCoilData.ini.

No Evaluation for FA1_FA2_FA3_FA4_FA5_FA6_FA7_FA8@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145262920

Series Loid for Evaluation 2:

Local 2.0.145262920

Number of images in series: 8

Evaluate image #1

Protocol 005_8Ch_FootAnkle_Sag_sn

Evaluation for FA1 specified in ServiceCoilData.ini.

Evaluation for coil element: FA1.

SN evaluation

New Evaluation mode: 0x743 specified

S/N Results

- Brightness = 820.2 [Specification from 500.0 to 2000.0]

- S/N Ratio = 109.0 [Specification > 80.0]

Applet Info:

Title: Series: 2, Image 1/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element FA1 are in Specification!

No Evaluation for FA1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_8Ch_FootAnkle_Sag_sn

Evaluation for FA2 specified in ServiceCoilData.ini.

Evaluation for coil element: FA2.

SN evaluation

New Evaluation mode: 0x534 specified

S/N Results

- Brightness = 1108.3 [Specification from 0.0 to 0.0]

- S/N Ratio = 103.6 [Specification > 100.0]

Applet Info:

Title: Series: 2, Image 2/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA2 are in Specification!
No Evaluation for FA2@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA3 specified in ServiceCoilData.ini.
Evaluation for coil element: FA3.
SN evaluation
New Evaluation mode: 0x743 specified
S/N Results
- Brightness = 1705.3 [Specification from 0.0 to 0.0]
- S/N Ratio = 172.8 [Specification > 100.0]
Applet Info:
Title: Series: 2, Image 3/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA3 are in Specification!
No Evaluation for FA3@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA4 specified in ServiceCoilData.ini.
Evaluation for coil element: FA4.
SN evaluation
New Evaluation mode: 0x745 specified
S/N Results
- Brightness = 1273.8 [Specification from 0.0 to 0.0]
- S/N Ratio = 121.9 [Specification > 100.0]
Applet Info:
Title: Series: 2, Image 4/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA4 are in Specification!
No Evaluation for FA4@1 specified in ServiceCoilData.ini.
Evaluate image #5
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA5 specified in ServiceCoilData.ini.
Evaluation for coil element: FA5.
SN evaluation
New Evaluation mode: 0x745 specified
S/N Results
- Brightness = 749.5 [Specification from 0.0 to 0.0]
- S/N Ratio = 101.0 [Specification > 80.0]
Applet Info:
Title: Series: 2, Image 5/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA5 are in Specification!
No Evaluation for FA5@1 specified in ServiceCoilData.ini.
Evaluate image #6
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA6 specified in ServiceCoilData.ini.
Evaluation for coil element: FA6.

SN evaluation
New Evaluation mode: 0x535 specified
S/N Results
- Brightness = 938.3 [Specification from 0.0 to 0.0]
- S/N Ratio = 109.7 [Specification > 90.0]
Applet Info:
Title: Series: 2, Image 6/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA6 are in Specification!
No Evaluation for FA6@1 specified in ServiceCoilData.ini.

Evaluate image #7
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA7 specified in ServiceCoilData.ini.
Evaluation for coil element: FA7.

SN evaluation
New Evaluation mode: 0x743 specified
S/N Results
- Brightness = 714.6 [Specification from 0.0 to 0.0]
- S/N Ratio = 217.4 [Specification > 150.0]
Applet Info:
Title: Series: 2, Image 7/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA7 are in Specification!
No Evaluation for FA7@1 specified in ServiceCoilData.ini.

Evaluate image #8
Protocol 005_8Ch_FootAnkle_Sag_sn
Evaluation for FA8 specified in ServiceCoilData.ini.
Evaluation for coil element: FA8.

SN evaluation
New Evaluation mode: 0x743 specified
S/N Results
- Brightness = 690.5 [Specification from 0.0 to 0.0]
- S/N Ratio = 260.1 [Specification > 150.0]
Applet Info:
Title: Series: 2, Image 8/8 of protocol 005_8Ch_FootAnkle_Sag_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element FA8 are in Specification!
No Evaluation for FA8@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished succesfully