Detected return code: 121

27641 QA

16Ch Body anterior Coil Check 2020-04-29 13:56:02

```
Start Coil Check
Performing Coil Check
Measure Protocol: 005 16Ch Body A BA1 BA2 Tra sn.edb.
Measured Series Loid #1:
Local 2.0.145240585
Measured Series Loid #2:
Local 2.0.145236212
Ordered Series Loid #1:
Local 2.0.145240585
Series Loid for Evaluation 1:
Local 2.0.145240585
Number of images in series: 2
Evaluate image #1
Protocol 005 16Ch_Body_A_BA1_BA2_Tra_sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A1A A1B A2B A2A A3A A3B A4B A4A specified in ServiceCoilData.ini.
Evaluation for coil element: AlA AlB A2B A2A A3A A3B A4B A4A.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness
                          = 1265.2 [Specification from 750.0 to 3000.0]
- S/N Ratio
                          = 66.9 [Specification >
Applet Info:
Title: Series: 1, Image 1/2 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A1A A1B A2B A2A A3A A3B A4B A4A are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A1A A1B A2B A2A A3A A3B A4B A4A@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: 32.000000 (specified), -32.000000 (image)
No Evaluation for A1A A1B A2B A2A A3A A3B A4B A4A specified in ServiceCoilData.ini.
Slice shift: 32.000000 (specified), -32.000000 (image)
No Evaluation for A1A A1B A2B A2A A3A A3B A4B A4A@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145236212
Series Loid for Evaluation 2:
Local 2.0.145236212
Number of images in series: 16
Evaluate image #1
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A1A specified in ServiceCoilData.ini.
Evaluation for coil element: A1A.
SN evaluation
```

```
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                       = 1525.7 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 91.5 [Specification >
Applet Info:
Title: Series: 2, Image 1/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element AlA are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A1A@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A1B specified in ServiceCoilData.ini.
Evaluation for coil element: A1B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
                         = 762.5 [Specification from 0.0 to 0.0]
- Brightness
- S/N Ratio
Applet Info:
                        = 66.9 [Specification > 50.0]
Title: Series: 2, Image 2/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A1B are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A1B@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A2B specified in ServiceCoilData.ini.
Evaluation for coil element: A2B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1430.8 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 84.7 [Specification > 50.0]
Applet Info:
Title: Series: 2, Image 3/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A2B are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A2B@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A2A specified in ServiceCoilData.ini.
Evaluation for coil element: A2A.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1335.4 [Specification from 0.0 to
- S/N Ratio
                        = 76.1 [Specification > 50.0]
Applet Info:
```

Title: Series: 2, Image 4/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb Data type: short. Width: 256 Height: 256 Results of coil-element A2A are in Specification! Slice shift: 32.000000 (specified), 32.000000 (image) No Evaluation for A2A@1 specified in ServiceCoilData.ini. Evaluate image #5 Protocol 005_16Ch_Body_A_BA1_BA2_Tra_sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A3A specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A3A@1 specified in ServiceCoilData.ini. Evaluate image #6 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A3B specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A3B@1 specified in ServiceCoilData.ini. Evaluate image #7 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A4B specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A4B@1 specified in ServiceCoilData.ini. Evaluate image #8 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A4A specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A4A@1 specified in ServiceCoilData.ini. Evaluate image #9 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A1A specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A1A@1 specified in ServiceCoilData.ini. Evaluate image #10 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A1B specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A1B@1 specified in ServiceCoilData.ini. Evaluate image #11 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A2B specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A2B@1 specified in ServiceCoilData.ini. Evaluate image #12 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A2A specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A2A@1 specified in ServiceCoilData.ini. Evaluate image #13 Protocol 005 16Ch Body A BA1 BA2 Tra sn Slice shift: -32.000000 (specified), -32.000000 (image) Evaluation for A3A specified in ServiceCoilData.ini. Evaluation for coil element: A3A. SN evaluation

```
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                       = 1656.8 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 102.2 [Specification >
Applet Info:
Title: Series: 2, Image 13/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A3A are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A3A@1 specified in ServiceCoilData.ini.
Evaluate image #14
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A3B specified in ServiceCoilData.ini.
Evaluation for coil element: A3B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
                         = 1484.8 [Specification from 0.0 to 0.0]
- Brightness
- S/N Ratio
Applet Info:
                        = 91.6 [Specification > 50.0]
Title: Series: 2, Image 14/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A3B are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A3B@1 specified in ServiceCoilData.ini.
Evaluate image #15
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A4B specified in ServiceCoilData.ini.
Evaluation for coil element: A4B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1298.4 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 92.3 [Specification > 50.0]
Applet Info:
Title: Series: 2, Image 15/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A4B are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A4B@1 specified in ServiceCoilData.ini.
Evaluate image #16
Protocol 005 16Ch Body A BA1 BA2 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A4A specified in ServiceCoilData.ini.
Evaluation for coil element: A4A.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1140.5 [Specification from 0.0 to
- S/N Ratio
                        = 76.4 [Specification > 50.0]
Applet Info:
```

```
Title: Series: 2, Image 16/16 of protocol 005 16Ch Body A BA1 BA2 Tra sn.edb
Data type:
             short
Width: 256
Height: 256
Results of coil-element A4A are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A4A@1 specified in ServiceCoilData.ini.
Detected return code: 122
Start Coil Check
Performing Coil Check
_____
Measure Protocol: 005 16Ch Body A_BA3_BA4_Tra_sn.edb.
Measured Series Loid #1:
Local 2.0.145242637
Measured Series Loid #2:
Local 2.0.145238266
Ordered Series Loid #1:
Local 2.0.145242637
Series Loid for Evaluation 1:
Local 2.0.145242637
Number of images in series: 2
Evaluate image #1
Protocol 005_16Ch_Body_A_BA3_BA4_Tra_sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A5A A5B A6B A6A A7A A7B A8B A8A specified in ServiceCoilData.ini.
Evaluation for coil element: A5A A5B A6B A6A A7A A7B A8B A8A.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness
                         = 1189.3 [Specification from 750.0 to 3000.0]
                         = 67.4 [Specification >
- S/N Ratio
                                                     50.01
Applet Info:
Title: Series: 1, Image 1/2 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A5A A5B A6B A6A A7A A7B A8B A8A are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A5A A5B A6B A6A A7A A7B A8B A8A@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: 32.000000 (specified), -32.000000 (image)
No Evaluation for A5A A5B A6B A6A A7A A7B A8B A8A specified in ServiceCoilData.ini.
Slice shift: 32.000000 (specified), -32.000000 (image)
No Evaluation for A5A A5B A6B A6A A7A A7B A8B A8A@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145238266
Series Loid for Evaluation 2:
Local 2.0.145238266
Number of images in series: 16
Evaluate image #1
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A5A specified in ServiceCoilData.ini.
Evaluation for coil element: A5A.
SN evaluation
```

```
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                       = 1443.7 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 105.2 [Specification >
Applet Info:
Title: Series: 2, Image 1/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A5A are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A5A@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A5B specified in ServiceCoilData.ini.
Evaluation for coil element: A5B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                         = 872.2 [Specification from 0.0 to 0.0]
- S/N Ratio
Applet Info:
                         = 70.7 [Specification > 50.0]
Title: Series: 2, Image 2/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A5B are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A5B@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A6B specified in ServiceCoilData.ini.
Evaluation for coil element: A6B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1250.4 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 81.5 [Specification > 50.0]
Applet Info:
Title: Series: 2, Image 3/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A6B are in Specification!
Slice shift: 32.000000 (specified), 32.000000 (image)
No Evaluation for A6B@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: 32.000000 (specified), 32.000000 (image)
Evaluation for A6A specified in ServiceCoilData.ini.
Evaluation for coil element: A6A.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1200.9 [Specification from 0.0 to
- S/N Ratio
                        = 93.8 [Specification > 50.0]
Applet Info:
```

Title: Series: 2, Image 4/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb Data type: short. Width: 256 Height: 256 Results of coil-element A6A are in Specification! Slice shift: 32.000000 (specified), 32.000000 (image) No Evaluation for A6A@1 specified in ServiceCoilData.ini. Evaluate image #5 Protocol 005_16Ch_Body_A_BA3_BA4_Tra_sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A7A specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A7A@1 specified in ServiceCoilData.ini. Evaluate image #6 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A7B specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A7B@1 specified in ServiceCoilData.ini. Evaluate image #7 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A8B specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A8B@1 specified in ServiceCoilData.ini. Evaluate image #8 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A8A specified in ServiceCoilData.ini. Slice shift: -32.000000 (specified), 32.000000 (image) No Evaluation for A8A@1 specified in ServiceCoilData.ini. Evaluate image #9 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A5A specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A5A@1 specified in ServiceCoilData.ini. Evaluate image #10 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A5B specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A5B@1 specified in ServiceCoilData.ini. Evaluate image #11 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A6B specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A6B@1 specified in ServiceCoilData.ini. Evaluate image #12 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A6A specified in ServiceCoilData.ini. Slice shift: 32.000000 (specified), -32.000000 (image) No Evaluation for A6A@1 specified in ServiceCoilData.ini. Evaluate image #13 Protocol 005 16Ch Body A BA3 BA4 Tra sn Slice shift: -32.000000 (specified), -32.000000 (image) Evaluation for A7A specified in ServiceCoilData.ini. Evaluation for coil element: A7A. SN evaluation

```
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                       = 1782.5 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 91.7 [Specification >
Applet Info:
Title: Series: 2, Image 13/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A7A are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A7A@1 specified in ServiceCoilData.ini.
Evaluate image #14
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A7B specified in ServiceCoilData.ini.
Evaluation for coil element: A7B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
                         = 1223.3 [Specification from 0.0 to 0.0]
- Brightness
- S/N Ratio
Applet Info:
                        = 83.1 [Specification > 50.0]
Title: Series: 2, Image 14/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A7B are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A7B@1 specified in ServiceCoilData.ini.
Evaluate image #15
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A8B specified in ServiceCoilData.ini.
Evaluation for coil element: A8B.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1028.2 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 74.4 [Specification > 50.0]
Applet Info:
Title: Series: 2, Image 15/16 of protocol 005 16Ch Body A BA3 BA4 Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element A8B are in Specification!
Slice shift: -32.000000 (specified), -32.000000 (image)
No Evaluation for A8B@1 specified in ServiceCoilData.ini.
Evaluate image #16
Protocol 005 16Ch Body A BA3 BA4 Tra sn
Slice shift: -32.000000 (specified), -32.000000 (image)
Evaluation for A8A specified in ServiceCoilData.ini.
Evaluation for coil element: A8A.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                        = 1133.5 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 82.2 [Specification > 50.0]
Applet Info:
```

Title: Series: 2, Image 16/16 of protocol 005_16Ch_Body_A_BA3_BA4_Tra_sn.edb

Data type: short

Width: 256 Height: 256

Results of coil-element A8A are in Specification! Slice shift: -32.000000 (specified), -32.000000 (image) No Evaluation for A8A@1 specified in ServiceCoilData.ini.

Results are in specification Program finished successfully