

## 27641 QA

### 16Ch Breast head Inv Coil Check 2020-04-29 14:35:01

---

Detected return code: 121

Start Coil Check

-----  
Performing Coil Check  
-----

Measure Protocol: 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb.

Measured Series Loid #1:

Local 2.0.145250857

Measured Series Loid #2:

Local 2.0.145246488

Ordered Series Loid #1:

Local 2.0.145250857

Series Loid for Evaluation 1:

Local 2.0.145250857

Number of images in series: 2

Evaluate image #1

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), 65.000000 (image)

Evaluation for BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8 specified in ServiceCoilData.ini.

Evaluation for coil element: BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1511.3 [Specification from 800.0 to 3200.0]

- S/N Ratio = 73.7 [Specification > 45.0]

Applet Info:

Title: Series: 1, Image 2/2 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)

No Evaluation for BL1\_BL2\_BL3\_BL4\_BL5\_BL6\_BL7\_BL8@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145246488

Series Loid for Evaluation 2:

Local 2.0.145246488

Number of images in series: 16

Evaluate image #1

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BL1 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BL1@1 specified in ServiceCoilData.ini.

Evaluate image #2  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: -40.000000 (specified), -40.000000 (image)  
Evaluation for BL2 specified in ServiceCoilData.ini.  
Evaluation for coil element: BL2.  
SN evaluation  
New Evaluation mode: 0x1634 specified  
S/N Results  
- Brightness = 422.0 [Specification from 0.0 to 0.0]  
- S/N Ratio = 35.9 [Specification > 22.0]  
Applet Info:  
Title: Series: 2, Image 2/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BL2 are in Specification!  
Slice shift: -40.000000 (specified), -40.000000 (image)  
No Evaluation for BL2@1 specified in ServiceCoilData.ini.  
Evaluate image #3  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL3 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL3@1 specified in ServiceCoilData.ini.  
Evaluate image #4  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL4 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL4@1 specified in ServiceCoilData.ini.  
Evaluate image #5  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL5 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL5@1 specified in ServiceCoilData.ini.  
Evaluate image #6  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL6 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL6@1 specified in ServiceCoilData.ini.  
Evaluate image #7  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL7 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL7@1 specified in ServiceCoilData.ini.  
Evaluate image #8  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL8 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BL8@1 specified in ServiceCoilData.ini.  
Evaluate image #9  
Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BL1 specified in ServiceCoilData.ini.  
Evaluation for coil element: BL1.  
SN evaluation  
New Evaluation mode: 0x1633 specified

## S/N Results

- Brightness = 259.7 [Specification from 0.0 to 0.0]  
- S/N Ratio = 16.1 [Specification > 5.0]

## Applet Info:

Title: Series: 2, Image 9/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

## Results of coil-element BL1 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BL1@1 specified in ServiceCoilData.ini.  
Evaluate image #10

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: -40.000000 (specified), 65.000000 (image)  
No Evaluation for BL2 specified in ServiceCoilData.ini.  
Slice shift: -40.000000 (specified), 65.000000 (image)  
No Evaluation for BL2@1 specified in ServiceCoilData.ini.  
Evaluate image #11

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BL3 specified in ServiceCoilData.ini.  
Evaluation for coil element: BL3.

## SN evaluation

New Evaluation mode: 0x1634 specified

## S/N Results

- Brightness = 1047.2 [Specification from 0.0 to 0.0]  
- S/N Ratio = 58.3 [Specification > 35.0]

## Applet Info:

Title: Series: 2, Image 11/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

## Results of coil-element BL3 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BL3@1 specified in ServiceCoilData.ini.  
Evaluate image #12

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BL4 specified in ServiceCoilData.ini.  
Evaluation for coil element: BL4.

## SN evaluation

New Evaluation mode: 0x1634 specified

## S/N Results

- Brightness = 868.6 [Specification from 0.0 to 0.0]  
- S/N Ratio = 60.6 [Specification > 30.0]

## Applet Info:

Title: Series: 2, Image 12/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

## Results of coil-element BL4 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BL4@1 specified in ServiceCoilData.ini.  
Evaluate image #13

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BL5 specified in ServiceCoilData.ini.  
Evaluation for coil element: BL5.  
SN evaluation

New Evaluation mode: 0x1634 specified

S/N Results

- Brightness = 571.3 [Specification from 0.0 to 0.0]  
- S/N Ratio = 41.2 [Specification > 30.0]

Applet Info:

Title: Series: 2, Image 13/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element BL5 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)

No Evaluation for BL5@1 specified in ServiceCoilData.ini.

Evaluate image #14

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), 65.000000 (image)

Evaluation for BL6 specified in ServiceCoilData.ini.

Evaluation for coil element: BL6.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 605.8 [Specification from 0.0 to 0.0]  
- S/N Ratio = 42.0 [Specification > 25.0]

Applet Info:

Title: Series: 2, Image 14/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element BL6 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)

No Evaluation for BL6@1 specified in ServiceCoilData.ini.

Evaluate image #15

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), 65.000000 (image)

Evaluation for BL7 specified in ServiceCoilData.ini.

Evaluation for coil element: BL7.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 907.3 [Specification from 0.0 to 0.0]  
- S/N Ratio = 48.2 [Specification > 30.0]

Applet Info:

Title: Series: 2, Image 15/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element BL7 are in Specification!

Slice shift: 65.000000 (specified), 65.000000 (image)

No Evaluation for BL7@1 specified in ServiceCoilData.ini.

Evaluate image #16

Protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), 65.000000 (image)

Evaluation for BL8 specified in ServiceCoilData.ini.

Evaluation for coil element: BL8.

SN evaluation

New Evaluation mode: 0x1634 specified

S/N Results

- Brightness = 809.4 [Specification from 0.0 to 0.0]  
- S/N Ratio = 55.3 [Specification > 30.0]

Applet Info:

Title: Series: 2, Image 16/16 of protocol 005\_16Ch\_Breast\_BL\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BL8 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BL8@1 specified in ServiceCoilData.ini.

-----  
Performing Coil Check  
-----

Measure Protocol: 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb.  
Measured Series Loid #1:  
Local 2.0.145252907  
Measured Series Loid #2:  
Local 2.0.145248542  
Ordered Series Loid #1:  
Local 2.0.145252907  
Series Loid for Evaluation 1:  
Local 2.0.145252907  
Number of images in series: 2  
Evaluate image #1  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8@1 specified in ServiceCoilData.ini.  
Evaluate image #2  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8.  
SN evaluation  
New Evaluation mode: 0x310 specified  
S/N Results  
- Brightness = 1492.2 [Specification from 800.0 to 3200.0]  
- S/N Ratio = 76.6 [Specification > 45.0]

Applet Info:  
Title: Series: 1, Image 2/2 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR1\_BR2\_BR3\_BR4\_BR5\_BR6\_BR7\_BR8@1 specified in ServiceCoilData.ini.  
Ordered Series Loid #2:  
Local 2.0.145248542  
Series Loid for Evaluation 2:  
Local 2.0.145248542  
Number of images in series: 16  
Evaluate image #1  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BR1 specified in ServiceCoilData.ini.  
Slice shift: 65.000000 (specified), -40.000000 (image)  
No Evaluation for BR1@1 specified in ServiceCoilData.ini.  
Evaluate image #2  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: -40.000000 (specified), -40.000000 (image)  
Evaluation for BR2 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR2.

SN evaluation

New Evaluation mode: 0x1634 specified

S/N Results

- Brightness = 362.1 [Specification from 0.0 to 0.0]  
- S/N Ratio = 31.9 [Specification > 22.0]

Applet Info:

Title: Series: 2, Image 2/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element BR2 are in Specification!

Slice shift: -40.000000 (specified), -40.000000 (image)

No Evaluation for BR2@1 specified in ServiceCoilData.ini.

Evaluate image #3

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR3 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR3@1 specified in ServiceCoilData.ini.

Evaluate image #4

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR4 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR4@1 specified in ServiceCoilData.ini.

Evaluate image #5

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR5 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR5@1 specified in ServiceCoilData.ini.

Evaluate image #6

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR6 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR6@1 specified in ServiceCoilData.ini.

Evaluate image #7

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR7 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR7@1 specified in ServiceCoilData.ini.

Evaluate image #8

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR8 specified in ServiceCoilData.ini.

Slice shift: 65.000000 (specified), -40.000000 (image)

No Evaluation for BR8@1 specified in ServiceCoilData.ini.

Evaluate image #9

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn

Slice shift: 65.000000 (specified), 65.000000 (image)

Evaluation for BR1 specified in ServiceCoilData.ini.

Evaluation for coil element: BR1.

SN evaluation

New Evaluation mode: 0x1633 specified

S/N Results

- Brightness = 191.6 [Specification from 0.0 to 0.0]  
- S/N Ratio = 11.3 [Specification > 5.0]

Applet Info:

Title: Series: 2, Image 9/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb

Data type: short  
Width: 256  
Height: 256

Results of coil-element BR1 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR1@1 specified in ServiceCoilData.ini.  
Evaluate image #10  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: -40.000000 (specified), 65.000000 (image)  
No Evaluation for BR2 specified in ServiceCoilData.ini.  
Slice shift: -40.000000 (specified), 65.000000 (image)  
No Evaluation for BR2@1 specified in ServiceCoilData.ini.  
Evaluate image #11  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR3 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR3.  
SN evaluation  
New Evaluation mode: 0x1634 specified  
S/N Results  
- Brightness = 1093.1 [Specification from 0.0 to 0.0]  
- S/N Ratio = 62.5 [Specification > 35.0]  
Applet Info:  
Title: Series: 2, Image 11/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR3 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR3@1 specified in ServiceCoilData.ini.  
Evaluate image #12  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR4 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR4.  
SN evaluation  
New Evaluation mode: 0x1634 specified  
S/N Results  
- Brightness = 560.8 [Specification from 0.0 to 0.0]  
- S/N Ratio = 40.7 [Specification > 30.0]  
Applet Info:  
Title: Series: 2, Image 12/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR4 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR4@1 specified in ServiceCoilData.ini.  
Evaluate image #13  
Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR5 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR5.  
SN evaluation  
New Evaluation mode: 0x1634 specified  
S/N Results  
- Brightness = 771.0 [Specification from 0.0 to 0.0]  
- S/N Ratio = 54.8 [Specification > 30.0]  
Applet Info:

Title: Series: 2, Image 13/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR5 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR5@1 specified in ServiceCoilData.ini.  
Evaluate image #14

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR6 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR6.  
SN evaluation  
New Evaluation mode: 0x310 specified

S/N Results  
- Brightness = 588.9 [Specification from 0.0 to 0.0]  
- S/N Ratio = 43.0 [Specification > 25.0]

Applet Info:

Title: Series: 2, Image 14/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR6 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR6@1 specified in ServiceCoilData.ini.  
Evaluate image #15

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR7 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR7.  
SN evaluation  
New Evaluation mode: 0x310 specified

S/N Results  
- Brightness = 914.9 [Specification from 0.0 to 0.0]  
- S/N Ratio = 51.2 [Specification > 30.0]

Applet Info:

Title: Series: 2, Image 15/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256

Results of coil-element BR7 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR7@1 specified in ServiceCoilData.ini.  
Evaluate image #16

Protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn  
Slice shift: 65.000000 (specified), 65.000000 (image)  
Evaluation for BR8 specified in ServiceCoilData.ini.  
Evaluation for coil element: BR8.  
SN evaluation  
New Evaluation mode: 0x1634 specified

S/N Results  
- Brightness = 761.9 [Specification from 0.0 to 0.0]  
- S/N Ratio = 51.8 [Specification > 30.0]

Applet Info:

Title: Series: 2, Image 16/16 of protocol 005\_16Ch\_Breast\_BR\_H\_Cor\_Inv\_sn.edb  
Data type: short  
Width: 256  
Height: 256



Results of coil-element BR8 are in Specification!  
Slice shift: 65.000000 (specified), 65.000000 (image)  
No Evaluation for BR8@1 specified in ServiceCoilData.ini.  
Results are in specification  
Program finished succesfully