

27641 QA

Spine Matrix Coil Check 2020-05-04 13:38:45

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_SpineMatrix_SP1_SP2_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145865959

Measured Series Loid #2:

Local 2.0.145848460

Ordered Series Loid #1:

Local 2.0.145865959

Series Loid for Evaluation 1:

Local 2.0.145865959

Number of images in series: 2

Evaluate image #1

Protocol 005_SpineMatrix_SP1_SP2_Tra_sn

No Evaluation for SP1_SP2 specified in ServiceCoilData.ini.

No Evaluation for SP1_SP2@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_SpineMatrix_SP1_SP2_Tra_sn

No Evaluation for SP1_SP2 specified in ServiceCoilData.ini.

No Evaluation for SP1_SP2@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145848460

Series Loid for Evaluation 2:

Local 2.0.145848460

Number of images in series: 4

Evaluate image #1

Protocol 005_SpineMatrix_SP1_SP2_Tra_sn

Slice shift: 70.000000 (specified), 70.000000 (image)

Evaluation for SP1 specified in ServiceCoilData.ini.

Evaluation for coil element: SP1.

SN evaluation

New Evaluation mode: 0x916 specified

S/N Results

- Brightness = 479.7 [Specification from 300.0 to 1200.0]

- S/N Ratio = 35.1 [Specification > 33.0]

Applet Info:

Title: Series: 2, Image 1/4 of protocol 005_SpineMatrix_SP1_SP2_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element SP1 are in Specification!

Slice shift: 70.000000 (specified), 70.000000 (image)

No Evaluation for SP1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_SpineMatrix_SP1_SP2_Tra_sn

Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP2 specified in ServiceCoilData.ini.
Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP2@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_SpineMatrix_SP1_SP2_Tra_sn
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP1 specified in ServiceCoilData.ini.
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP1@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_SpineMatrix_SP1_SP2_Tra_sn
Slice shift: -70.000000 (specified), -70.000000 (image)
Evaluation for SP2 specified in ServiceCoilData.ini.
Evaluation for coil element: SP2.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results
- Brightness = 514.1 [Specification from 0.0 to 0.0]
- S/N Ratio = 38.0 [Specification > 33.0]
Applet Info:
Title: Series: 2, Image 4/4 of protocol 005_SpineMatrix_SP1_SP2_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element SP2 are in Specification!
Slice shift: -70.000000 (specified), -70.000000 (image)
No Evaluation for SP2@1 specified in ServiceCoilData.ini.
Detected return code: 122

Start Coil Check

Performing Coil Check

Measure Protocol: 005_SpineMatrix_SP3_SP4_Tra_sn.edb.
Measured Series Loid #1:
Local 2.0.145869037
Measured Series Loid #2:
Local 2.0.145851534
Ordered Series Loid #1:
Local 2.0.145869037
Series Loid for Evaluation 1:
Local 2.0.145869037
Number of images in series: 2
Evaluate image #1
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
No Evaluation for SP3_SP4 specified in ServiceCoilData.ini.
No Evaluation for SP3_SP4@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
No Evaluation for SP3_SP4 specified in ServiceCoilData.ini.
No Evaluation for SP3_SP4@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145851534
Series Loid for Evaluation 2:
Local 2.0.145851534
Number of images in series: 4

Evaluate image #1
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
Slice shift: 70.000000 (specified), 70.000000 (image)
Evaluation for SP3 specified in ServiceCoilData.ini.
Evaluation for coil element: SP3.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results
- Brightness = 498.1 [Specification from 0.0 to 0.0]
- S/N Ratio = 35.4 [Specification > 33.0]
Applet Info:
Title: Series: 2, Image 1/4 of protocol 005_SpineMatrix_SP3_SP4_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element SP3 are in Specification!
Slice shift: 70.000000 (specified), 70.000000 (image)
No Evaluation for SP3@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP4 specified in ServiceCoilData.ini.
Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP4@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP3 specified in ServiceCoilData.ini.
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP3@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_SpineMatrix_SP3_SP4_Tra_sn
Slice shift: -70.000000 (specified), -70.000000 (image)
Evaluation for SP4 specified in ServiceCoilData.ini.
Evaluation for coil element: SP4.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results
- Brightness = 531.5 [Specification from 0.0 to 0.0]
- S/N Ratio = 38.6 [Specification > 33.0]
Applet Info:
Title: Series: 2, Image 4/4 of protocol 005_SpineMatrix_SP3_SP4_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element SP4 are in Specification!
Slice shift: -70.000000 (specified), -70.000000 (image)
No Evaluation for SP4@1 specified in ServiceCoilData.ini.
Detected return code: 123

Start Coil Check

Performing Coil Check

Measure Protocol: 005_SpineMatrix_SP5_SP6_Tra_sn.edb.
Measured Series Loid #1:

Local 2.0.145871089
Measured Series Loid #2:
Local 2.0.145854608
Ordered Series Loid #1:
Local 2.0.145871089
Series Loid for Evaluation 1:
Local 2.0.145871089
Number of images in series: 2
Evaluate image #1
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
No Evaluation for SP5_SP6 specified in ServiceCoilData.ini.
No Evaluation for SP5_SP6@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
No Evaluation for SP5_SP6 specified in ServiceCoilData.ini.
No Evaluation for SP5_SP6@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145854608
Series Loid for Evaluation 2:
Local 2.0.145854608
Number of images in series: 4
Evaluate image #1
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
Slice shift: 70.000000 (specified), 70.000000 (image)
Evaluation for SP5 specified in ServiceCoilData.ini.
Evaluation for coil element: SP5.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results
- Brightness = 512.1 [Specification from 0.0 to 0.0]
- S/N Ratio = 35.0 [Specification > 33.0]
Applet Info:
Title: Series: 2, Image 1/4 of protocol 005_SpineMatrix_SP5_SP6_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element SP5 are in Specification!
Slice shift: 70.000000 (specified), 70.000000 (image)
No Evaluation for SP5@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP6 specified in ServiceCoilData.ini.
Slice shift: -70.000000 (specified), 70.000000 (image)
No Evaluation for SP6@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP5 specified in ServiceCoilData.ini.
Slice shift: 70.000000 (specified), -70.000000 (image)
No Evaluation for SP5@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_SpineMatrix_SP5_SP6_Tra_sn
Slice shift: -70.000000 (specified), -70.000000 (image)
Evaluation for SP6 specified in ServiceCoilData.ini.
Evaluation for coil element: SP6.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results

- Brightness = 528.7 [Specification from 0.0 to 0.0]
- S/N Ratio = 36.3 [Specification > 33.0]

Applet Info:

Title: Series: 2, Image 4/4 of protocol 005_SpineMatrix_SP5_SP6_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element SP6 are in Specification!

Slice shift: -70.000000 (specified), -70.000000 (image)

No Evaluation for SP6@1 specified in ServiceCoilData.ini.

Detected return code: 124

Start Coil Check

Performing Coil Check

Measure Protocol: 005_SpineMatrix_SP7_SP8_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145873141

Measured Series Loid #2:

Local 2.0.145867934

Ordered Series Loid #1:

Local 2.0.145873141

Series Loid for Evaluation 1:

Local 2.0.145873141

Number of images in series: 2

Evaluate image #1

Protocol 005_SpineMatrix_SP7_SP8_Tra_sn

No Evaluation for SP7_SP8 specified in ServiceCoilData.ini.

No Evaluation for SP7_SP8@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_SpineMatrix_SP7_SP8_Tra_sn

No Evaluation for SP7_SP8 specified in ServiceCoilData.ini.

No Evaluation for SP7_SP8@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145867934

Series Loid for Evaluation 2:

Local 2.0.145867934

Number of images in series: 4

Evaluate image #1

Protocol 005_SpineMatrix_SP7_SP8_Tra_sn

Slice shift: 2.000000 (specified), -120.000000 (image)

No Evaluation for SP7 specified in ServiceCoilData.ini.

Slice shift: 2.000000 (specified), -120.000000 (image)

No Evaluation for SP7@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_SpineMatrix_SP7_SP8_Tra_sn

Slice shift: -120.000000 (specified), -120.000000 (image)

Evaluation for SP8 specified in ServiceCoilData.ini.

Evaluation for coil element: SP8.

SN evaluation

New Evaluation mode: 0x916 specified

S/N Results

- Brightness = 529.7 [Specification from 0.0 to 0.0]

- S/N Ratio = 39.1 [Specification > 33.0]

Applet Info:

Title: Series: 2, Image 2/4 of protocol 005_SpineMatrix_SP7_SP8_Tra_sn.edb

Data type: short
Width: 256
Height: 256

Results of coil-element SP8 are in Specification!
Slice shift: -120.000000 (specified), -120.000000 (image)
No Evaluation for SP8@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_SpineMatrix_SP7_SP8_Tra_sn
Slice shift: 2.000000 (specified), 2.000000 (image)
Evaluation for SP7 specified in ServiceCoilData.ini.
Evaluation for coil element: SP7.
SN evaluation
New Evaluation mode: 0x916 specified
S/N Results
- Brightness = 559.0 [Specification from 0.0 to 0.0]
- S/N Ratio = 37.3 [Specification > 33.0]
Applet Info:
Title: Series: 2, Image 3/4 of protocol 005_SpineMatrix_SP7_SP8_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element SP7 are in Specification!
Slice shift: 2.000000 (specified), 2.000000 (image)
No Evaluation for SP7@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_SpineMatrix_SP7_SP8_Tra_sn
Slice shift: -120.000000 (specified), 2.000000 (image)
No Evaluation for SP8 specified in ServiceCoilData.ini.
Slice shift: -120.000000 (specified), 2.000000 (image)
No Evaluation for SP8@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished succesfully