

27641 QA

PAMatrix feet Coil Check 2020-04-29 13:11:23

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_PAMatrix_PL1_PR1_F_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145223975

Measured Series Loid #2:

Local 2.0.145219792

Ordered Series Loid #1:

Local 2.0.145223975

Series Loid for Evaluation 1:

Local 2.0.145223975

Number of images in series: 2

Evaluate image #1

Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn

No Evaluation for PL1_PR1 specified in ServiceCoilData.ini.

No Evaluation for PL1_PR1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn

No Evaluation for PL1_PR1 specified in ServiceCoilData.ini.

No Evaluation for PL1_PR1@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145219792

Series Loid for Evaluation 2:

Local 2.0.145219792

Number of images in series: 4

Evaluate image #1

Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn

Slice shift: 10.000000 (specified), -10.000000 (image)

No Evaluation for PL1 specified in ServiceCoilData.ini.

Slice shift: 10.000000 (specified), -10.000000 (image)

No Evaluation for PL1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn

Slice shift: -10.000000 (specified), -10.000000 (image)

Evaluation for PR1 specified in ServiceCoilData.ini.

Evaluation for coil element: PR1.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1474.4 [Specification from 750.0 to 3000.0]

- S/N Ratio = 80.1 [Specification > 52.0]

Applet Info:

Title: Series: 2, Image 2/4 of protocol 005_PAMatrix_PL1_PR1_F_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element PR1 are in Specification!
Slice shift: -10.000000 (specified), -10.000000 (image)
No Evaluation for PR1@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn
Slice shift: 10.000000 (specified), 10.000000 (image)
Evaluation for PL1 specified in ServiceCoilData.ini.
Evaluation for coil element: PL1.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness = 1312.3 [Specification from 500.0 to 3000.0]
- S/N Ratio = 72.9 [Specification > 52.0]
Applet Info:
Title: Series: 2, Image 3/4 of protocol 005_PAMatrix_PL1_PR1_F_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element PL1 are in Specification!
Slice shift: 10.000000 (specified), 10.000000 (image)
No Evaluation for PL1@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_PAMatrix_PL1_PR1_F_Tra_sn
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR1 specified in ServiceCoilData.ini.
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR1@1 specified in ServiceCoilData.ini.
Detected return code: 122

Start Coil Check

Performing Coil Check

Measure Protocol: 005_PAMatrix_PL2_PR2_F_Tra_sn.edb.
Measured Series Loid #1:
Local 2.0.145226027
Measured Series Loid #2:
Local 2.0.145221844
Ordered Series Loid #1:
Local 2.0.145226027
Series Loid for Evaluation 1:
Local 2.0.145226027
Number of images in series: 2
Evaluate image #1
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
No Evaluation for PL2_PR2 specified in ServiceCoilData.ini.
No Evaluation for PL2_PR2@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
No Evaluation for PL2_PR2 specified in ServiceCoilData.ini.
No Evaluation for PL2_PR2@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145221844
Series Loid for Evaluation 2:
Local 2.0.145221844
Number of images in series: 4

Evaluate image #1
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL2 specified in ServiceCoilData.ini.
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL2@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
Slice shift: -10.000000 (specified), -10.000000 (image)
Evaluation for PR2 specified in ServiceCoilData.ini.
Evaluation for coil element: PR2.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness = 1583.1 [Specification from 0.0 to 0.0]
- S/N Ratio = 85.1 [Specification > 57.0]
Applet Info:
Title: Series: 2, Image 2/4 of protocol 005_PAMatrix_PL2_PR2_F_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element PR2 are in Specification!
Slice shift: -10.000000 (specified), -10.000000 (image)
No Evaluation for PR2@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
Slice shift: 10.000000 (specified), 10.000000 (image)
Evaluation for PL2 specified in ServiceCoilData.ini.
Evaluation for coil element: PL2.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness = 1366.9 [Specification from 0.0 to 0.0]
- S/N Ratio = 85.1 [Specification > 57.0]
Applet Info:
Title: Series: 2, Image 3/4 of protocol 005_PAMatrix_PL2_PR2_F_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element PL2 are in Specification!
Slice shift: 10.000000 (specified), 10.000000 (image)
No Evaluation for PL2@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_PAMatrix_PL2_PR2_F_Tra_sn
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR2 specified in ServiceCoilData.ini.
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR2@1 specified in ServiceCoilData.ini.
Detected return code: 123

Start Coil Check

Performing Coil Check

Measure Protocol: 005_PAMatrix_PL3_PR3_F_Tra_sn.edb.
Measured Series Loid #1:

```

Local 2.0.145228079
Measured Series Loid #2:
Local 2.0.145223896
Ordered Series Loid #1:
Local 2.0.145228079
Series Loid for Evaluation 1:
Local 2.0.145228079
Number of images in series: 2
Evaluate image #1
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
No Evaluation for PL3_PR3 specified in ServiceCoilData.ini.
No Evaluation for PL3_PR3@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
No Evaluation for PL3_PR3 specified in ServiceCoilData.ini.
No Evaluation for PL3_PR3@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145223896
Series Loid for Evaluation 2:
Local 2.0.145223896
Number of images in series: 4
Evaluate image #1
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL3 specified in ServiceCoilData.ini.
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL3@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
Slice shift: -10.000000 (specified), -10.000000 (image)
Evaluation for PR3 specified in ServiceCoilData.ini.
Evaluation for coil element: PR3.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness                = 1095.8 [Specification from 0.0 to 0.0]
- S/N Ratio                 = 111.2 [Specification > 75.0]
Applet Info:
Title: Series: 2, Image 2/4 of protocol 005_PAMatrix_PL3_PR3_F_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element PR3 are in Specification!
Slice shift: -10.000000 (specified), -10.000000 (image)
No Evaluation for PR3@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
Slice shift: 10.000000 (specified), 10.000000 (image)
Evaluation for PL3 specified in ServiceCoilData.ini.
Evaluation for coil element: PL3.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness                = 1092.7 [Specification from 0.0 to 0.0]
- S/N Ratio                 = 100.5 [Specification > 75.0]
Applet Info:
Title: Series: 2, Image 3/4 of protocol 005_PAMatrix_PL3_PR3_F_Tra_sn.edb
Data type: short
Width: 256

```

Height: 256

Results of coil-element PL3 are in Specification!
Slice shift: 10.000000 (specified), 10.000000 (image)
No Evaluation for PL3@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_PAMatrix_PL3_PR3_F_Tra_sn
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR3 specified in ServiceCoilData.ini.
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR3@1 specified in ServiceCoilData.ini.
Detected return code: 124

Start Coil Check

Performing Coil Check

Measure Protocol: 005_PAMatrix_PL4_PR4_F_Tra_sn.edb.
Measured Series Loid #1:
Local 2.0.145230131
Measured Series Loid #2:
Local 2.0.145225950
Ordered Series Loid #1:
Local 2.0.145230131
Series Loid for Evaluation 1:
Local 2.0.145230131
Number of images in series: 2
Evaluate image #1
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
No Evaluation for PL4_PR4 specified in ServiceCoilData.ini.
No Evaluation for PL4_PR4@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
No Evaluation for PL4_PR4 specified in ServiceCoilData.ini.
No Evaluation for PL4_PR4@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145225950
Series Loid for Evaluation 2:
Local 2.0.145225950
Number of images in series: 4
Evaluate image #1
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL4 specified in ServiceCoilData.ini.
Slice shift: 10.000000 (specified), -10.000000 (image)
No Evaluation for PL4@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
Slice shift: -10.000000 (specified), -10.000000 (image)
Evaluation for PR4 specified in ServiceCoilData.ini.
Evaluation for coil element: PR4.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness = 840.7 [Specification from 0.0 to 0.0]
- S/N Ratio = 92.1 [Specification > 65.0]
Applet Info:
Title: Series: 2, Image 2/4 of protocol 005_PAMatrix_PL4_PR4_F_Tra_sn.edb

Data type: short
Width: 256
Height: 256

Results of coil-element PR4 are in Specification!
Slice shift: -10.000000 (specified), -10.000000 (image)
No Evaluation for PR4@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
Slice shift: 10.000000 (specified), 10.000000 (image)
Evaluation for PL4 specified in ServiceCoilData.ini.
Evaluation for coil element: PL4.
SN evaluation
New Evaluation mode: 0x310 specified
S/N Results
- Brightness = 753.5 [Specification from 0.0 to 0.0]
- S/N Ratio = 90.8 [Specification > 65.0]
Applet Info:
Title: Series: 2, Image 3/4 of protocol 005_PAMatrix_PL4_PR4_F_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element PL4 are in Specification!
Slice shift: 10.000000 (specified), 10.000000 (image)
No Evaluation for PL4@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005_PAMatrix_PL4_PR4_F_Tra_sn
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR4 specified in ServiceCoilData.ini.
Slice shift: -10.000000 (specified), 10.000000 (image)
No Evaluation for PR4@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished succesfully