

27641 QA

Head Matrix Coil Check 2020-04-29 16:23:58

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_HeadMatrix_HE1_HE2_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145278374

Measured Series Loid #2:

Local 2.0.145274218

Ordered Series Loid #1:

Local 2.0.145278374

Series Loid for Evaluation 1:

Local 2.0.145278374

Number of images in series: 1

Evaluate image #1

Protocol 005_HeadMatrix_HE1_HE2_Tra_sn

No Evaluation for HE1_HE2 specified in ServiceCoilData.ini.

No Evaluation for HE1_HE2@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145274218

Series Loid for Evaluation 2:

Local 2.0.145274218

Number of images in series: 2

Evaluate image #1

Protocol 005_HeadMatrix_HE1_HE2_Tra_sn

Evaluation for HE1 specified in ServiceCoilData.ini.

Evaluation for coil element: HE1.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 951.9 [Specification from 400.0 to 3000.0]

- S/N Ratio = 81.0 [Specification > 75.0]

Applet Info:

Title: Series: 2, Image 1/2 of protocol 005_HeadMatrix_HE1_HE2_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element HE1 are in Specification!

No Evaluation for HE1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_HeadMatrix_HE1_HE2_Tra_sn

Evaluation for HE2 specified in ServiceCoilData.ini.

Evaluation for coil element: HE2.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1350.8 [Specification from 0.0 to 0.0]

- S/N Ratio = 98.8 [Specification > 80.0]

Applet Info:

Title: Series: 2, Image 2/2 of protocol 005_HeadMatrix_HE1_HE2_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element HE2 are in Specification!

No Evaluation for HE2@1 specified in ServiceCoilData.ini.

Performing Coil Check

Measure Protocol: 005_HeadMatrix_HE3_HE4_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145280424

Measured Series Loid #2:

Local 2.0.145276272

Ordered Series Loid #1:

Local 2.0.145280424

Series Loid for Evaluation 1:

Local 2.0.145280424

Number of images in series: 1

Evaluate image #1

Protocol 005_HeadMatrix_HE3_HE4_Tra_sn

No Evaluation for HE3_HE4 specified in ServiceCoilData.ini.

No Evaluation for HE3_HE4@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145276272

Series Loid for Evaluation 2:

Local 2.0.145276272

Number of images in series: 2

Evaluate image #1

Protocol 005_HeadMatrix_HE3_HE4_Tra_sn

Evaluation for HE3 specified in ServiceCoilData.ini.

Evaluation for coil element: HE3.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 884.4 [Specification from 0.0 to 0.0]

- S/N Ratio = 65.6 [Specification > 50.0]

Applet Info:

Title: Series: 2, Image 1/2 of protocol 005_HeadMatrix_HE3_HE4_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element HE3 are in Specification!

No Evaluation for HE3@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_HeadMatrix_HE3_HE4_Tra_sn

Evaluation for HE4 specified in ServiceCoilData.ini.

Evaluation for coil element: HE4.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1048.9 [Specification from 0.0 to 0.0]

- S/N Ratio = 71.2 [Specification > 55.0]

Applet Info:

Title: Series: 2, Image 2/2 of protocol 005_HeadMatrix_HE3_HE4_Tra_sn.edb

Data type: short

Width: 256
Height: 256

Results of coil-element HE4 are in Specification!
No Evaluation for HE4@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished succesfully