

## Quality Assurance (04-May-2020, 10:51:32, #141814)

---

**Coil Check Spine\_32 Success (04-May-2020, 10:51:33 » 04-May-2020, 11:11:10)**

---

### Comments (0)

---

#	Timestamp	Comment
No comments for this function.		

### Summary

---

Summary

	Status
098_Spine_32_SP1_Tra	OK
098_Spine_32_SP2_Tra	OK
098_Spine_32_SP3_Tra	OK
098_Spine_32_SP4_Tra	OK
098_Spine_32_SP5_Tra	OK
098_Spine_32_SP6_Tra	OK
098_Spine_32_SP7_Tra	OK
098_Spine_32_SP8_Tra	OK
Program Result	Success

### Details

---

#### Protocol: 098\_Spine\_32\_SP1\_Tra

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

#### Coil: S11

---

Signal



S/N

Parabola Fit	
Value	Min
41.8	29.0

Brightness

Parabola Fit		
Value	Min	Max
363	180	700


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	304	V

**Coil: S12**

---

Signal



S/N

Parabola Fit	
Value	Min

35.2	26.0
------	------

Brightness

Parabola Fit
360


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	304	V

**Coil: S13**

---

Signal



S/N

Parabola Fit	
Value	Min
36.8	26.0

Brightness

Parabola Fit
325


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	304	V

**Coil: S14**

---

Signal



S/N

Parabola Fit	
Value	Min
51.6	38.0

Brightness

Parabola Fit
541

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	304	V

**Protocol: 098\_Spine\_32\_SP2\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S21**

---

Signal



S/N

Parabola Fit	
Value	Min
41.4	30.0

Brightness

Parabola Fit
462


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S22**

---

Signal



S/N

Parabola Fit	
Value	Min
32.2	25.0

Brightness

Parabola Fit
290


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S23**

---

Signal



S/N

Parabola Fit	
Value	Min
32.2	23.0

Brightness

Parabola Fit
326


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S24**

---

Signal



S/N

Parabola Fit	
Value	Min
41.8	30.0

Brightness

Parabola Fit
353

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Protocol: 098\_Spine\_32\_SP3\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S31**

---

Signal



S/N

Parabola Fit	
Value	Min
38.6	29.0

Brightness

Parabola Fit
366


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S32**

---

Signal



S/N

Parabola Fit	
Value	Min
36.5	26.0



Brightness

Parabola Fit
353


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S33**

---

Signal



S/N

Parabola Fit	
Value	Min
36.8	26.0

Brightness

Parabola Fit
348


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S34**

---

Signal



S/N

Parabola Fit	
Value	Min
49.1	38.0

Brightness

Parabola Fit
519

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Protocol: 098\_Spine\_32\_SP4\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S41**

---

Signal



S/N

Parabola Fit	
Value	Min
39.9	30.0

Brightness

Parabola Fit
480


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S42**

---

Signal



S/N

Parabola Fit	
Value	Min
32.6	25.0

Brightness

Parabola Fit
301


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S43**

---

Signal



S/N

Parabola Fit	
Value	Min
32.8	23.0

Brightness

Parabola Fit
328


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S44**

---

Signal



S/N

Parabola Fit	
Value	Min
41.0	30.0

Brightness

Parabola Fit
357

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Protocol: 098\_Spine\_32\_SP5\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S51**

---

Signal



S/N

Parabola Fit	
Value	Min
36.2	29.0

Brightness

Parabola Fit
367


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S52**

---

Signal



S/N

Parabola Fit	
Value	Min
35.5	26.0

Brightness

Parabola Fit
369


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S53**

---

Signal



S/N

Parabola Fit	
Value	Min
34.6	26.0

Brightness

Parabola Fit
343


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Coil: S54**

---

Signal



S/N

Parabola Fit	
Value	Min
49.0	38.0

Brightness

Parabola Fit
546

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	302	V

**Protocol: 098\_Spine\_32\_SP6\_Tra**

---

Coil Info


Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S61**

---

Signal





S/N

Parabola Fit	
Value	Min
40.1	30.0

Brightness

Parabola Fit
461


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	301	V

**Coil: S62**

---

Signal



S/N

Parabola Fit	
Value	Min
31.8	25.0

Brightness

Parabola Fit
297


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	301	V

**Coil: S63**

---

Signal



S/N

Parabola Fit	
Value	Min
32.2	23.0

Brightness

Parabola Fit
333


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	301	V

**Coil: S64**

---

Signal



S/N

Parabola Fit	
Value	Min
40.8	30.0

Brightness

Parabola Fit
368

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	301	V

**Protocol: 098\_Spine\_32\_SP7\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S71**

---

Signal



S/N

Parabola Fit	
Value	Min
37.9	29.0

Brightness

Parabola Fit
356


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S72**

---

Signal



S/N

Parabola Fit	
Value	Min
34.2	26.0

Brightness

Parabola Fit
353


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S73**

---

Signal



S/N

Parabola Fit	
Value	Min
35.5	26.0

Brightness

Parabola Fit
339


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S74**

---

Signal



S/N

Parabola Fit	
Value	Min
51.9	38.0

Brightness

Parabola Fit
541

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Protocol: 098\_Spine\_32\_SP8\_Tra**

---


Coil Info

Coil ID	Serial Number	Part Number
Spine_32	30999	10496544

**Coil: S81**

---

Signal



S/N

Parabola Fit	
Value	Min
24.0	17.0

Brightness

Parabola Fit
215


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S82**

---

Signal



S/N

Parabola Fit	
Value	Min
24.3	14.0

Brightness

Parabola Fit
167


Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S83**

---

Signal



S/N

Parabola Fit	
Value	Min
25.6	14.0

Brightness

Parabola Fit
200

Homogeneity and Transmitter Amplitude


Variable	Value	Unit
Transmitter Amplitude	303	V

**Coil: S84**

---

Signal





S/N

Parabola Fit	
Value	Min
30.3	20.0

Brightness

Parabola Fit
195

Homogeneity and Transmitter Amplitude

Variable	Value	Unit
Transmitter Amplitude	303	V