

27641 QA

TxRx CP Extremity Coil Check 2020-04-29 15:54:02

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_Extremity_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145270106

No valid multi seriesloid.

Measured Series Loid #1:

Local 2.0.145270106

Ordered Series Loid #1:

Local 2.0.145270106

Series Loid for Evaluation 1:

Local 2.0.145270106

Number of images in series: 1

Evaluate image #1

Protocol 005_Extremity_Tra_sn

Evaluation for EX specified in ServiceCoilData.ini.

Evaluation for coil element: EX.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1363.7 [Specification from 700.0 to 2800.0]

- S/N Ratio = 86.9 [Specification > 80.0]

Applet Info:

Title: Series: 1, Image 1/1 of protocol 005_Extremity_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element EX are in Specification!

No Evaluation for EX@1 specified in ServiceCoilData.ini.

Results are in specification

Program finished succesfully