

27641 QA

Neck Matrix Coil Check 2020-04-29 16:10:06

Detected return code: 121

Start Coil Check

Performing Coil Check

Measure Protocol: 005_NeckMatrix_all_Tra_sn.edb.

Measured Series Loid #1:

Local 2.0.145276320

Measured Series Loid #2:

Local 2.0.145272164

Ordered Series Loid #1:

Local 2.0.145276320

Series Loid for Evaluation 1:

Local 2.0.145276320

Number of images in series: 1

Evaluate image #1

Protocol 005_NeckMatrix_all_Tra_sn

No Evaluation for NE1_NE2 specified in ServiceCoilData.ini.

No Evaluation for NE1_NE2@1 specified in ServiceCoilData.ini.

Ordered Series Loid #2:

Local 2.0.145272164

Series Loid for Evaluation 2:

Local 2.0.145272164

Number of images in series: 2

Evaluate image #1

Protocol 005_NeckMatrix_all_Tra_sn

Evaluation for NE1 specified in ServiceCoilData.ini.

Evaluation for coil element: NE1.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 1057.3 [Specification from 500.0 to 2000.0]

- S/N Ratio = 70.5 [Specification > 60.0]

Applet Info:

Title: Series: 2, Image 1/2 of protocol 005_NeckMatrix_all_Tra_sn.edb

Data type: short

Width: 256

Height: 256

Results of coil-element NE1 are in Specification!

No Evaluation for NE1@1 specified in ServiceCoilData.ini.

Evaluate image #2

Protocol 005_NeckMatrix_all_Tra_sn

Evaluation for NE2 specified in ServiceCoilData.ini.

Evaluation for coil element: NE2.

SN evaluation

New Evaluation mode: 0x310 specified

S/N Results

- Brightness = 907.7 [Specification from 0.0 to 0.0]

- S/N Ratio = 78.1 [Specification > 60.0]
Applet Info:
Title: Series: 2, Image 2/2 of protocol 005_NeckMatrix_all_Tra_sn.edb
Data type: short
Width: 256
Height: 256

Results of coil-element NE2 are in Specification!
No Evaluation for NE2@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished succesfully