## 27641 QA

## 8Ch Wrist Coil Check 2020-04-29 15:40:21

```
Detected return code: 121
Start Coil Check
Performing Coil Check
-----
Measure Protocol: 005 8Ch Wrist Tra sn.edb.
Measured Series Loid #1:
Local 2.0.145270249
Measured Series Loid #2:
Local 2.0.145265998
Ordered Series Loid #1:
Local 2.0.145270249
Series Loid for Evaluation 1:
Local 2.0.145270249
Number of images in series: 1
Evaluate image #1
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR1 WR2 WR3 WR4 WR5 WR6 WR7 WR8 specified in ServiceCoilData.ini.
Evaluation for coil element: WR1 WR2 WR3 WR4 WR5 WR6 WR7 WR8.
SN evaluation
New Evaluation mode: 0x340 specified
S/N Results
- Brightness
                         = 1078.5 [Specification from 400.0 to 2000.0]
                         = 167.4 [Specification > 120.0]
- S/N Ratio
Applet Info:
Title: Series: 1, Image 1/1 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR1 WR2 WR3 WR4 WR5 WR6 WR7 WR8 are in Specification!
No Evaluation for WR1 WR2 WR3 WR4 WR5 WR6 WR7 WR8@1 specified in ServiceCoilData.ini.
Ordered Series Loid #2:
Local 2.0.145265998
Series Loid for Evaluation 2:
Local 2.0.145265998
Number of images in series: 8
Evaluate image #1
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR1 specified in ServiceCoilData.ini.
Evaluation for coil element: WR1.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                         = 533.6 [Specification from 0.0 to
                                                                  0.01
- S/N Ratio
                         = 88.3 [Specification > 60.0]
Applet Info:
Title: Series: 2, Image 1/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type:
            short
Width: 256
Height: 256
```

```
Results of coil-element WR1 are in Specification!
No Evaluation for WR1@1 specified in ServiceCoilData.ini.
Evaluate image #2
Protocol 005_8Ch_Wrist_Tra_sn
Evaluation for WR2 specified in ServiceCoilData.ini.
Evaluation for coil element: WR2.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
- Brightness
                         = 442.7 [Specification from 0.0 to 0.0]
- S/N Ratio
                         = 76.0 [Specification > 50.0]
Applet Info:
Title: Series: 2, Image 2/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR2 are in Specification!
No Evaluation for WR2@1 specified in ServiceCoilData.ini.
Evaluate image #3
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR3 specified in ServiceCoilData.ini.
Evaluation for coil element: WR3.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
- Brightness
                         = 438.5 [Specification from 0.0 to 0.0]
- S/N Ratio
                         = 75.4 [Specification >
Applet Info:
Title: Series: 2, Image 3/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR3 are in Specification!
No Evaluation for WR3@1 specified in ServiceCoilData.ini.
Evaluate image #4
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR4 specified in ServiceCoilData.ini.
Evaluation for coil element: WR4.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
                        = 632.8 [Specification from 0.0 to 0.0]
- Brightness
                         = 100.4 [Specification > 65.0]
- S/N Ratio
Applet Info:
Title: Series: 2, Image 4/8 of protocol 005 8Ch Wrist Tra sn.edb
             short
Data type:
Width: 256
Height: 256
Results of coil-element WR4 are in Specification!
No Evaluation for WR4@1 specified in ServiceCoilData.ini.
Evaluate image #5
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR5 specified in ServiceCoilData.ini.
Evaluation for coil element: WR5.
SN evaluation
New Evaluation mode: 0x1643 specified
S/N Results
- Brightness
                         = 580.2 [Specification from 0.0 to 0.0]
```

```
- S/N Ratio
                       = 105.1 [Specification > 65.0]
Applet Info:
Title: Series: 2, Image 5/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR5 are in Specification!
No Evaluation for WR5@1 specified in ServiceCoilData.ini.
Evaluate image #6
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR6 specified in ServiceCoilData.ini.
Evaluation for coil element: WR6.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
- Brightness
                        = 573.3 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 102.8 [Specification > 65.0]
Applet Info:
Title: Series: 2, Image 6/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR6 are in Specification!
No Evaluation for WR601 specified in ServiceCoilData.ini.
Evaluate image #7
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR7 specified in ServiceCoilData.ini.
Evaluation for coil element: WR7.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
- Brightness
                         = 432.2 [Specification from 0.0 to 0.0]
                         = 73.6 [Specification > 40.0]
- S/N Ratio
Applet Info:
Title: Series: 2, Image 7/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR7 are in Specification!
No Evaluation for WR7@1 specified in ServiceCoilData.ini.
Evaluate image #8
Protocol 005 8Ch Wrist Tra sn
Evaluation for WR8 specified in ServiceCoilData.ini.
Evaluation for coil element: WR8.
SN evaluation
New Evaluation mode: 0x1642 specified
S/N Results
- Brightness
                         = 477.6 [Specification from 0.0 to 0.0]
- S/N Ratio
                        = 81.6 [Specification > 55.0]
Title: Series: 2, Image 8/8 of protocol 005 8Ch Wrist Tra sn.edb
Data type: short
Width: 256
Height: 256
Results of coil-element WR8 are in Specification!
No Evaluation for WR8@1 specified in ServiceCoilData.ini.
Results are in specification
Program finished successfully
```